

Orgo-Life® Rife Frequency Generator UL-9051



Brand: Orgo-Life® Technology

Product Code: UL-9051

Availability: 7

Weight: 0.00kg

Dimensions: 0.00cm x 0.00cm x 0.00cm

Price: \$189.95

Ex Tax: \$189.95

Short Description

Orgo-Life® Rife Frequency Generator UL-9051—Simple , fast & Ideal for beginners doing Bio-Resonance

The **Orgo-Life® UL-9051 Frequency Generator** is an enhanced, stabilized version designed specifically for the Orgo-Life® ecosystem. It works in full harmony with **all vibrational units, quantum modules, coils, antennas, energetic plates, and wellness test systems.**

Description

UL-9051 Rife System – 1 Hz to 500 kHz Frequency Generator + UL7013 Interface + Full Orgo-Life® Accessories

quantum modules, coils, antennas, energetic plates, and wellness test systems.

Everything is **100% non-medical**, intended **exclusively for wellness purposes**, fully compliant with **Law 31**.

This system allows you to perform:

- personal tests
- harmonization frequency sessions
- wellness-exploration tests

- vibrational experiences to restore inner balance
- reduction of “bad vibes”
- calming internal tensions
- regaining energetic stability in the body

The generator is officially delivered with **all Orgo-Life® cables and interfaces**, ready to connect to every UL device, spiral, module, plate, coil and Orgo-Life® accessory.

1 Official Orgo-Life® Kit Content

Each system includes:

— Orgo-Life® UL-7013 Interface

Allows connection to:

- UL modules
- spirals
- headsets
- plates
- discs
- antennas
- copper tubes
- Orgo-Life® quantum accessories

— Mini 1/8 Interface (included)

A compact universal adapter allowing connection to **all Orgo-Life® units**, including older generations.

— Orgo-Life® Tesla Coil Antenna (your choice)

- 50 mm
 - 100 mm
 - 150 mm

Each coil includes:

- ? dedicated cable
- ? Orgo-Life® adapter
- ? UL-7013 interface connection support

— Orgo-Life® Copper Tubes (2 units)

- Length: 6 inches
 - Diameter: $\frac{3}{4}$ inch
 - Material: 100% pure copper
- Used for hand-held frequency tests or harmonization sessions.

— Orgo-Life® BNC Alligator Clip Cables

Used to connect:

- UL modules
- plates
- copper tubes
- antennas
- various vibrational test setups

— Additional Coil Cable (50 / 100 / 150 mm depending on your choice)

Designed to transmit frequencies with stability.

— High-stability USB power cable (3C certified)

2 General Description of the UL-9051 Orgo-Life® Generator

This DDS (Direct Digital Synthesis) generator can produce **over 15 vibrational waveform types**, all compatible with Orgo-Life® systems.

Available Waveforms:

- Sine
 - Square
 - Triangle
 - Sawtooth + Reverse Sawtooth
 - ECG
 - Schumann (7.83 Hz default)
 - Gaussian
 - Lorentz
 - Sinc
 - Exponential rise
 - Exponential decay
 - Half-wave
 - Full-wave
 - Noise
- ...and all standard international waveform variations used in vibrational

testing.

Technical Specifications:

- **Frequency range:**
 - Sine: 0.01 Hz ? 500 kHz
 - Other waveforms: 0.01 Hz ? 100 kHz
- **Resolution:** 0.01 Hz (ultra-precise)
- **Max amplitude:** $\pm 10\text{V}$ (Peak-to-Peak)
- **Impedance:** 50 ohms
- **DC Offset:** $\pm 10\text{V}$ (switchable)
- **Power supply:** USB DC 5V

Integrated Filters

Designed to produce cleaner signals, especially for:

- high-frequency sine waves
- pulsed signals
- Rodin coils / Orgo-Life® antennas

Typical Orgo-Life® Applications:

- Harmonization
 - Vibrational cleansing
 - Personal quantum exploration
 - Emotional stabilization (wellness only, non-medical)
 - Subtle meridian realignment
 - Mental calming
 - Exploratory frequency testing

3 Controls and Functions

MODE — Changes the waveform

CURSOR — Selects the digit to edit

+ / ? — Adjusts the frequency

RUN/STOP — Starts or stops output

DCOFFSET — Activates/deactivates DC component

FILTER — Cleans the output signal

AMP — Adjusts signal amplitude

4 Steps for Using With Orgo-Life®

1. Plug the USB power cable into the generator.
2. Connect the BNC cable to the UL-7013 interface.
3. Connect your modules:
 - coil antenna
 - pure copper tubes ($\frac{3}{4}$ inch $\times$ 6 inches)
 - Tesla spiral 50 mm, 100 mm or 150 mm
 - quantum plate
4. Select your waveform.
5. Enter a frequency (e.g., 7.83 Hz, 432 Hz, 963 Hz).
6. Press **RUN**.

5 Wellness Application (Law 31 Compliant)

This generator is used for:

- centering oneself
- reducing negative energy or internal discomfort
- regaining a stable state
- supporting deep relaxation
- improving vibrational sensitivity
- exploring frequencies that help “clear the mind”
- harmonizing intentions and emotions

Always for **wellness purposes only**, never medical, in accordance with **Law 31**.

6 Mandatory Legal Notice (Law 31)

Orgo-Life® devices are intended exclusively for wellness, harmonization and energetic relaxation.

They **do not diagnose, treat, or prevent any disease**.

They do not replace any healthcare professional.

For any physical or psychological condition, consult a licensed medical specialist.

TECHNICAL SPECIFICATIONS

MODEL: UL-9051

DESCRIPTION:

- 3 $\times$ downloadable PDF booklets (on Orgo-Life.com) with product description (FR/EN)
- 1 $\times$ USB power block
- 1 $\times$ USB charging cable
- 1 $\times$ alligator cable (optional) with 2 $\times$ $\frac{3}{4}$ inch pure copper pipes

- 1 × UL-7013 interface cable to plug Orgo-Life® units
- 1 × Tesla antenna 50 mm, 100 mm or 150 mm with cable and adapter plug
- 1 × BNC to mini 1/8 cable for FI headphone (very comfortable)

VOLTAGE:

110 volts or 220 volts – all devices are universal.

All devices come with 220 V adapters (for countries outside North America), compatible with French/European current, and work on 240 V as well as 110 V.

PACKAGING:

35 cm × 25 cm × 15 cm – approx. 3 kg

BARCODE:

787416647851

INSTRUCTIONS:

Comes with French & English instruction manuals in PDF (downloadable on Orgo-Life.com)

DISCLAIMER:

Our devices do not claim to cure any disease or illness.

In the event of a serious medical condition, please seek the help and advice of a qualified medical professional.

The potential health benefits have not been evaluated by national medical licensing bodies.

MANUFACTURER:

Original Health Distribution Ltd., Canada

Email: original.health.net@gmail.com


PRODUCT SAFETY:

Instructions are downloadable on the Orgo-Life.com website (PDF document).

PRODUCTION PARTNERS:

OrgoLifeTechnology, with support from Renewable Free Energy Foundation, Magog, QC, Canada.

Product Gallery



Orgo-Life®
3 years Limited Warranty
Garantie Limitée de 3 ans

Product Code :
UL-8051

QUANTUM ENERGY MEDBED TECHNOLOGY®

Ideal for Alternative Practitioner
Idéal pour les praticiens alternatifs

Orgo-Life.com & Med-bed.co

Orgo-Life®
3 years Limited Warranty
Garantie Limitée de 3 ans

Product Code :
UL-8051

Use the copper pipes options to do special Rife frequency work on your body

Ideal for Alternative Practitioner
Idéal pour les praticiens alternatifs

Orgo-Life.com & Med-bed.co

Orgo-Life®
3 years Limited Warranty
Garantie Limitée de 3 ans

Product Code :
UL-8051

Option without copper pipes & its cable
With Tesla Antenna 50 mm + cables

Ideal for Alternative Practitioner
Idéal pour les praticiens alternatifs

Orgo-Life.com & Med-bed.co

Orgo-Life®
3 years Limited Warranty
Garantie Limitée de 3 ans

Product Code :
UL-8051

-Ideal for Rife long distance energy healing applications, two coils sit over your head, to the coronal chakra (center of receptance, that sets the Rife frequency on the left example with the Tesla Antenna 5, 10 or 15 cm greater energy output)

Ideal for Alternative Practitioner
Idéal pour les praticiens alternatifs

Orgo-Life.com & Med-bed.co

Orgo-Life®
3 years Limited Warranty
Garantie Limitée de 3 ans

Product Code :
UL-8051

Compatible with all Orgo-life® Units
Also compatible with Quantum Boosters
UL-8055 (1000%) & UL-8050 (300%)

Ideal for Alternative Practitioner
Idéal pour les praticiens alternatifs

Orgo-Life.com & Med-bed.co

Orgo-Life®
3 years Limited Warranty
Garantie Limitée de 3 ans

Product Code :
UL-8051

Option without copper pipes & its cable
With Tesla Antenna 100 mm + cables

Ideal for Alternative Practitioner
Idéal pour les praticiens alternatifs

Orgo-Life.com & Med-bed.co



